


# analyser

the art of measuring



User Manual  
FlowMeter F1, F2

IMT.Analytics

IMT Analytics AG  
Gewerbstrasse 8  
9470 Buchs (SG)  
Switzerland

[www.imtanalytics.com](http://www.imtanalytics.com)

# Table of contents

1	Introduction	5
2	Intended use	6
3	Safety instructions	7
3.1	Representation of hazards, cautions and notes	7
3.2	Personnel	7
3.3	Responsibility and guarantee	7
3.4	Service life	8
4	Specifications	9
4.1	Measurement parameters	9
4.2	Principle of flow measurement	10
4.3	Special functions	10
4.4	Communication interfaces	10
4.5	Physical data	10
4.6	Calibration by user	11
4.7	Operating data	11
4.8	Gas standards for flow and volume measurement	12
4.9	Power supply	12
4.10	Battery operation	13
4.11	Directives and approvals	13
4.12	Device labels and symbols	14
4.13	Feature Set	15
5	Start-up	16
5.1	Individual parts in the packaging	16
5.2	Power supply	17
5.3	Mechanical connectors	17
5.4	Electrical interfaces	20
6	Operation	22
6.1	Switching the device on and off	22
6.2	Acoustic feedback	22
6.3	User control gestures	22
6.4	Screens	24
6.5	Menu	28
7	Measuring volume	41
7.1	General principle	41
7.2	Finding the correct trigger values	41
8	Servicing and care	42
8.1	Guidelines for servicing and care	42
8.2	Notes about changing parts	42
8.3	Preventive cleaning and servicing routines	42
8.4	Cleaning	43
8.5	Trouble shooting	43
8.6	Contact	45

9	Accessories and spare parts	46
9.1	Ordering address	46
9.2	Parts	46
10	Disposal	48
11	Measurement parameters and units	49
11.1	Measurement parameters and units	49
12	Wireless communication	50
13	Appendix	52
13.1	Abbreviations and glossary	52
13.2	Index	54

# 1 Introduction

## Validity

This documentation is valid for the product: FlowMeter F1 and F2

You will find the name FlowMeter on the type label on the back of your device.

## Software and firmware version

This documentation applies to the following versions:

FlowMeter firmware – Version 1.2.000 or higher

In the case of older or newer versions there may be minor discrepancies in relation to this User Manual.

## Designations used in this User Manual

Buttons such as **Power** and indicators on the display such as **Settings** are printed in boldface italics.

## References to pages and sections

For references to pages and references, e.g. ([→4.5 Physical data](#)), the symbol ([→XY](#)) is used.

## Version information

Issue date of this User Manual: **Release 05, 2023-08**

Subject to technical modifications without notice.

## 2 Intended use

This product is intended for testing and calibration purposes of medical – or other devices or systems that generate Air / O<sub>2</sub> / CO<sub>2</sub><sup>1</sup> / N<sub>2</sub> gas flows and / or gas pressures.

The device is intended for users who are trained in the use of measuring devices and who can carry out tests, repairs, calibrations, maintenance and service work on medical or other devices. No dedicated training on this specific device is required.

The device can be used in:

- Hospital service departments
- Clinics
- Medical and non-medical devices manufacturing facilities
- Workshop / testing / laboratory facilities of independent service companies performing service and repair for medical and non-medical devices

FlowMeter is intended for use in an indoor test laboratory environment. It may only be used outside the nursing sector. It must not be used directly on patients or devices that are connected to patients. The FlowMeter is intended for over-the-counter sale.

FlowMeter is the solution for measurements in the following areas:

- Flow ( $\pm 300$  L/min)
- Volume (0 L – 500 000 L)
- Differential pressure ( $\pm 250$  mbar)
- Absolute pressure in the Flow channel (0.5 – 2 bar)
- Temperature in the Flow channel (-10 – 50 °C)
- Humidity in the Flow channel (0 – 100 % RH, non-condensing)



FlowMeter is a measuring instrument for testing and calibrating purpose on medical and non-medical devices or systems. It must not be used for patient monitoring. FlowMeter must not be used directly on patients or devices that are connected to patients.

This product is intended to be used at elevations of up to 5000 m operation altitude in buildings.

<sup>1</sup> Only available for the model FlowMeter F2

## 3 Safety instructions

### 3.1 Representation of hazards, cautions and notes

This User Manual uses the representation below to specifically draw attention to residual risks during intended use and emphasize important technical requirements.



Information and/or instructions and prohibitions to prevent damage of any kind.



Use CO<sub>2</sub> and other non-breathable gases only in well ventilated rooms. Increased ambient concentration of non-breathable gases can lead to unconsciousness or death.

### 3.2 Personnel



Work on and with FlowMeter may only be performed by persons who have undergone appropriate general technical training and have the necessary experience.

### 3.3 Responsibility and guarantee

The manufacturer accepts no responsibility or guarantee and will be exempt from liability claims accordingly if the operator or any third parties:

- fail to use the device in accordance with its intended use
- disregard the specifications
- tamper with the device in any way (conversions, modifications, etc.)
- operate the device with accessories that are not listed in the associated sets of product documentation



- Report any serious incident that has occurred in relation to the medical device to the manufacturer (→8.6.2 [Technical support](#)) and the authority having jurisdiction in their locale
- The device has user replaceable parts (→8.3 [Preventive cleaning and servicing routines](#))
- Negligence of warning and error messages can lead to damage of the device
- The device including the measurement channel shall not be exposed to high level of volatile organic compounds (VOC). Doing so may lead to permanent offset of the humidity sensor.
- Only use original packaging in good condition for transport or storage
- The device must not be placed in a pressure chamber
- The device is not intended to measure flammable or combustible gases
- The device is not intended to measure liquids
- The device must not be connected to a telecom network
- The device must only be connect to approved equipment with CE, CSA, UL or other equivalent safety mark to ensure double insulation



Although the device meets high quality and safety standards and has been constructed and tested according to the current state of the art, it is not possible to rule out the risk of injuries with serious consequences if the device is not used in compliance with the intended use (improperly) or is misused. Therefore, please read this User Manual carefully and keep this documentation in a readily accessible place close to your device.

### 3.4 Service life

The typical service life of the device is specified as 10 (ten) years, provided it is handled properly in accordance with this User Manual.



## 4 Specifications

### 4.1 Measurement parameters

Parameter		Value
Flow (Air / O <sub>2</sub> / N <sub>2</sub> ) <sup>1</sup> (→4.13 Feature Set)	Range Accuracy	±300 L/min FlowMeter F1: ±1.9 % or ±0.05 L/min FlowMeter F2: ±1.7 % or ±0.05 L/min
Flow (CO <sub>2</sub> ) <sup>2</sup>	Range Accuracy	±100 L/min FlowMeter F2: ±3.5 % or ±0.05 L/min
Volume	Range Accuracy	0 L – 500 000 L FlowMeter F1: ± 2.25 % or ± 1mL (flow > 2.5 L/min) FlowMeter F2: ± 2 % or ± 1mL (flow > 2.5 L/min)
Differential pressure <sup>2</sup>	Range Accuracy	±250 mbar ±0.5 % or ±0.15 mbar
Absolute pressure in the Flow channel	Range Accuracy	0.5 – 2 bar ±10 mbar
Temperature <sup>3</sup>	Range Accuracy	-10 – 50 °C ±1 °C
Humidity <sup>2</sup>	Range Accuracy	0 – 100 % RH* < 20 % RH: ±5 % RH* 20 % – 80 % RH: ±3 % RH* > 80 % RH: ±5 % RH*

#### Accuracy valid

- between 10 °C and 30 °C (channel temperature)
- between 950 hPa and 1250 hPa
- zero offset calibration performed after warm-up
- for measurements performed within 4 h after last zero offset calibration at same ambient conditions

<sup>1</sup> Flow accuracy specified is valid under the above mentioned general conditions and for a flow in positive flow direction.

- Add 0.05 % of reading per 1 °C outside the range of 10 °C to 30 °C (channel temperature)
- Add 0.005 % of reading per 10 hPa above 1250 hPa
- Add 0.01 % of reading per 10 hPa below 950 hPa
- Add 0.05 % of reading per 10 % oxygen concentration above 21 %
- Add 0.2 % of reading for negative flow direction

Flow accuracy stated with gas (channel) temperature and ambient temperature within 10 °C of each other. Add 0.05 % of reading per 1 °C above 10 °C temperature difference.

<sup>2</sup> Only available for the model FlowMeter F2

<sup>3</sup> At flow rate ≥ 50 L/min, screen brightness ≤ 30 %, battery not charging

\* Non-condensing




All flow and volume specifications are based on STP gas standard at 21.1 °C and 1013.25 mbar.

For all specifications with absolute and relative tolerances, the greater value applies (for example: ±1.9 % of reading or ±0.05 L/min absolute, whichever is greater).

**4.2 Principle of flow measurement**

Flow in the flow channel is determined by differential pressure measurement. To build up differential pressure  $\Delta p$ , a measuring disc is used to provide flow resistance. For flow rates  $> 50 \text{ L/min}$  the differential pressure  $\Delta p_2$  is used. For lower flows, the Flow-Meter uses  $\Delta p_1$  to improve the accuracy at lower flow rates.



$$\Delta p = c_1 \cdot \eta \cdot Q + c_2 \cdot \rho \cdot Q^2$$

$\eta$ : dynamic viscosity of the gas [Pa\*s]  
 $\rho$ : gas density [kg/m<sup>3</sup>]  
 $c_1, c_2$ : device-specific constants (channel geometry)

**Dynamic viscosity**

The viscosity of a medium is its resistance to the flow and shear of the current. Viscosity is extremely dependent on temperature. The viscosity of a medium is slightly dependent on the pressure and moisture content of the medium.

**Density**

Density is the unit for the mass per unit volume of the medium. Density is very dependent on pressure and temperature.

The impact of ambient conditions is hence the reason why flow is occasionally converted to standard conditions.

(→4.8 Gas standards for flow and volume measurement)

**4.3 Special functions**

Automatic battery operation in the event of a power failure.

**4.4 Communication interfaces**

- USB-A: for data recording and software updates
- USB-C: no communication, used only for charging
- RJ-10 (serial port) to stream measurement values and set up the device

**4.5 Physical data**

Weight: 350 g  
 Size (l x w x h): 200 x 80 x 60 mm  
 Gas types: Air, O<sub>2</sub>, CO<sub>2</sub><sup>1</sup>, N<sub>2</sub> and mixtures: Air/O<sub>2</sub>

<sup>1</sup> Only available for the model FlowMeter F2

**4.6 Calibration by user**

Offset calibration of the pressure and flow sensors (→[6.5.8 Zero offset calibration](#)).

**4.7 Operating data**

Degree of protection, against water and dust, according to IEC 60529: IP 20

**The two devices FlowMeter F1 and FlowMeter F2 have 3 different modes:**

- Battery operated
- Mains supplied / battery charging
- Switched off / storage

Please note the different ambient condition for the different modes. The allowed conditions for the Flow channel and the device enclosure are different.

**Battery operated**

Environment conditions	Allowed range
Temperature (flow channel / enclosure)	-10 – 50 °C (14 – 122 °F)
Air humidity (enclosure)	10 % – 90 % RH*
Air humidity (flow channel)	0 % – 100 % RH*
Absolute pressure in flow channel	50 kPa – 200 kPa
Atmospheric pressure	54 kPa – 120 kPa

**Mains operated**

Environment conditions	Allowed range
Temperature (flow channel / enclosure)	5 – 40 °C (41 – 104 °F)
Air humidity (enclosure)	10 % – 90 % RH*
Air humidity (flow channel)	0 % – 100 % RH*
Absolute pressure in flow channel	50 kPa – 200 kPa
Atmospheric pressure	54 kPa – 120 kPa

**Switched off / storage**

Environment conditions	Allowed range
Temperature (flow channel / enclosure)	-10 – 60 °C (14 – 140 °F)
Air humidity (flow channel / enclosure)	5 – 95 % RH*
Absolute pressure in flow channel	50 kPa – 200 kPa
Atmospheric pressure	54 kPa – 110 kPa

\* Non-condensing

**4.8 Gas standards for flow and volume measurement**

FlowMeter converts the flow and volume readings measured in the device to match the conditions of the standard selected. The following gas standards are supported by FlowMeter:

Gas standard		Temperature	Pressure	Relative humidity
Ambient Temperature and Pressure	ATP	Current gas temperature	Current ambient pressure	Current gas humidity
Channel Temperature and Pressure	CTP	Current gas temperature	Current channel pressure	Current gas humidity
Standard Conditions USA	STP	21.1 °C (70 °F)	1013.25 mbar (760 mmHg)	0 %
Body Temperature and (Ambient) Pressure Saturated according to ISO 80601-2-12:2011	BTPS-A	37 °C (99 °F)	Current ambient pressure	100 %



Please refer to [→11 Measurement parameters and units](#). Here, you will also find the conversion factors for the units of measurement.

**4.9 Power supply**

Input voltage of the power supply: 100 – 240 VAC (±10 %), 50 – 60 Hz, 0.6 A  
 Output voltage of the power supply: 5 VDC, 3.0 A  
 FlowMeter input voltage: 5 VDC, 2.5 A

The original power supply from IMT Analytics must be used for trouble-free and reliable operation. It is power tested and approved according to IEC 62368-1.

#### 4.10 Battery operation

Operating time during battery operation: at least 4 hours.

##### Charging the battery

A complete charging process takes about 4 hours in normal conditions. The service life of the battery is extended if the battery is not completely discharged. Do not store the device with a discharged battery. Optimal state of charge for long term storage is 80 %.



The battery is not user replaceable and must not be replaced with a different type. Doing so may lead to an explosion, burn or fire hazard.



The device indicates visually and audibly when the battery must be charged. Please do not store the device if the battery is in a depleted state.

**Caution: Complete depletion can damage the battery beyond repair!**

4

#### 4.11 Directives and approvals

- IEC 61010-1
- IEC 61326-1
- CAN/CSA-C22.2 No. 61010-1-12
- UL 61010-1 3<sup>rd</sup> Edition



The device is classified as Pollution Degree 2.

The device is classified as Overvoltage Category II.

##### For the USB connection

V1.1 is used (12 Mbit)

##### For RJ-10 connection

If the device is to be actuated via the RS-232 interface, your dealer will be pleased to provide you with a detailed protocol.



The device is not intended for use outside a building.






















##### Simplified EU Declaration of Conformity







Hereby, IMT Analytics AG declares that the radio equipment types FlowMeter F1 and FlowMeter F2 are in compliance with Directive 2014/53/EU.

The full text of the EU declaration of conformity is available at the following internet address: <https://www.imtanalytics.com/Extranet/Folders>

4.12 Device labels and symbols

The following labels and symbols can be found on FlowMeter, packaging or accessories:

Symbol	Description
	Warning notice
	Power button
	Alternating current
	Class II equipment
	For indoor use only
	Direct current
	Power input
	USB Port to connect an USB flash drive
	Serial interface
	Conforms to CE Directives and Regulations
	Safety Mark for North America for FlowMeter
	Mark for proof of product compliance to North American safety standards
	Manufacturer
	Date of Manufacture
	The operating instructions should be considered when operating the device
	Serial Number
	Waste Electrical and Electronic Equipment
	Conformity with Low Voltage Directive 2014/35/EU
	Mark for United Kingdom Conformity Assessment
	This equipment contains specified radio equipment that has been certified to the Technical Regulation Conformity Certification under the Radio Law
 UDI Data  Matrix Code 	Unique Device Identification (01) GTIN-No. (11) Production Date (21) Serial No.

Symbol	Description
	Catalog number
	Keep dry
	Protect from sunlight
	Temperature for transport and storage
	Humidity range for transport and storage
	Atmospheric pressure for transport and storage

4

### 4.13 Feature Set

The two devices FlowMeter F1 and FlowMeter F2 differ from each other in terms of the available features.

The models can be easily distinguished by the color of the top cover.

	FlowMeter F1	FlowMeter F2
Top cover	Black	Blue
Flow measurement	±300 L/min ±1.9 % or ±0.05 L/min	±300 L/min ±1.7 % or ±0.05 L/min
CO <sub>2</sub> flow measurement	not available	available
Volume measurement	available	available
Absolute pressure in flow channel	available	available
Differential pressure ports	not available	available
Temperature measurement	available	available
Humidity measurement	not available	available
Bluetooth	available	available
Statistic screen	not available	available

## 5 Start-up

### 5.1 Individual parts in the packaging

Picture	Description
 A blue and white rectangular flow meter device with a digital display and two ports on the front.	FlowMeter
 A white power supply unit with a USB cable and four different types of power adapters (two-prong, three-prong, and two-prong with ground).	Power supply
 A clear, cylindrical plastic protection filter with a wider base and a narrower top section.	Protection Filter RT019
 A stack of documents, including a black folder and several white certificates with blue text and logos.	Calibration certificate



## 5.2 Power supply

The power supply socket is located at the rear of FlowMeter.



The device can be disconnected from the mains by disconnecting the power supply. The power supply should therefore be easily accessible.

### 5.2.1 Supply voltage

The mains voltage of the power supply is 100 – 240 VAC at 50 – 60 Hz.



Before switching on, make sure the operating voltage of the power supply is compatible with the local mains voltage. You will find this information on the rating plate on the back of the power supply.

## 5.3 Mechanical connectors

### 5.3.1 Protection Filter RT019


To protect the device against contaminants and particles in the flow channel, the Protection Filter RT019, must be used for flow measurements. Connect the Protection Filter RT019 to the flow channel inlet port before connecting the device under test.



Particles of dirt in the air can clog the measuring system and thus lead to incorrect measurements. The Protection Filter RT019 must be checked regularly (→8.3 Preventive cleaning and servicing routines).

### 5.3.2 Flow channel

The flow channel can be used to perform gas flow measurements from -300 L/min to 300 L/min. The channel includes pressure, temperature and humidity<sup>1</sup> sensors.



Measuring range & accuracy details: see chapter ([→4.1 Measurement parameters](#))



When working with relatively high humidity, ensure that condensation does not form in the unit! Water can irreparably damage the sensors!




The device including the measurement channel shall not be exposed to high level of volatile organic compounds (VOC). Doing so may lead to permanent offset of the humidity sensor.

<sup>1</sup> Only available for the model FlowMeter F2

### 5.3.3 Differential pressure<sup>1</sup>

The differential pressure connectors can be used for differential pressure measurements.



Differential pressure ports

Measuring range:  $\pm 250$  mbar

Accuracy:  $\pm 0.5$  % or  $\pm 0.15$  mbar

<sup>1</sup> Only available for the model FlowMeter F2

## 5.4 Electrical interfaces

### 5.4.1 USB connection for charging

The USB-C port is used to connect the power supply. The port is located at the rear of the device.



Rating: 5 V, 2.5 A

The device must only be connect to approved equipment with CE, CSA, UL or other equivalent safety mark to ensure double insulation.

### 5.4.2 RS-232

The RS-232 interface is used for remote monitoring or control.




Actuation of the RS-232 port takes place via a special RS-232 cable. If the device shall be actuated via the RS-232 interface, your dealer will be pleased to provide you with a detailed protocol.

FlowMeter pin assignment (RJ-10 connector):

Pin 1	NC
Pin 2	RxD (Input of FlowMeter)
Pin 3	TxD (Output of FlowMeter)
Pin 4	GND

### 5.4.3 USB-A

The USB-A host connector can be used to connect a USB flash drive for data recording and software update.



Rating: 5 V, 0.5 A

The USB flash drive must comply with the following requirements.

- File system: FAT32
- Volume size: 0.5 – 32 GByte

USB hubs are not supported.

## 6 Operation

### 6.1 Switching the device on and off

The device is switched on and off using the power button on the rear of the device.






The device provides a beep as acoustic feedback as following:


Event	Duration of beep	Number of beeps
Power on/off	Short	1
Battery state of charge falls below 20 %	Short	1
Battery state of charge falls below 10 %	Short	2
Battery state of charge falls below 5 %	Short	3
Device shuts down due to empty battery	Long	1
Info message	Short	1
Warning message	Medium	1
Error message	Medium	3

### 6.3 User control gestures

To operate the touch screen optimally and easily, use the gestures listed here:

Gesture		Action	Function
	Tap	Tap an element once	<ul style="list-style-type: none"> <li>Selection of menu items</li> <li>Change between Numerical view and Chart view</li> </ul>
	Long press (> 0.5 sec.)	Tap and hold for one second	<ul style="list-style-type: none"> <li>Editing readings and curves</li> </ul>
	Swipe	Drag your finger across the screen	<ul style="list-style-type: none"> <li>Change views on the measurement screen</li> </ul>

If FlowMeter is switched on, the Measurement screen appears.



Two screens can be selected by swiping your finger left or right:


No.	Description
1	Measurement screen
2	Statistic screen

Note: Statistic Screen is only available on FlowMeter F2

6.4 Screens



6.4.1 Measurement screen

The measurement screen is the starting point to operate the FlowMeter. The software version described in this User Manual displays the menu items listed below:



No.	Description	Action
1	Value tiles Three different tiles can be configured individually with the following parameters: <ul style="list-style-type: none"> <li>Measurement value (Updates every 1 sec.)</li> <li>Unit</li> </ul> <div style="display: flex; justify-content: space-around; font-size: small;"> <div style="text-align: center;">Humidity <b>55.0</b> %</div> <div style="text-align: center;">Temperature <b>21.4</b> °C</div> <div style="text-align: center;">Volume <b>535.0</b> ml</div> </div> <p>If a value is outside of the specified measurement range, the background color will change to red:</p> <div style="text-align: center; background-color: #f08080; padding: 5px; font-size: small;">                     Temperature <b>59.0</b> °C                 </div>	Long press: The parameter configuration menu opens
2	Chart tile parameters Shows the properties of the parameters: <ul style="list-style-type: none"> <li>Name</li> <li>Unit</li> </ul>	Long press: The parameter configuration menu opens
3	Menu button Shows the Menu	Tap: Overview of all sub menus opens (→6.5 Menu)
4	Gas standard Shows current configured gas standard	Tap: Selection of all gas standards opens
5	Gas type Shows current configured gas type. Only for Air/O <sub>2</sub> : selection for O <sub>2</sub> concentration is available.	Tap: Selection of all gas types opens
6	Recording status Indicates if recording is active: <ul style="list-style-type: none"> <li>Shown: Recording active</li> <li>Hidden: Recording inactive</li> </ul>	Tap: Recording menu opens
7	Bluetooth status Indicates the status of the bluetooth connection	




No.	Description	Action
8	<p>Battery level</p> <p>Indicates the battery level:</p> <ul style="list-style-type: none"> <li>• Empty battery: No bar</li> <li>• Full battery: Full bar</li> <li>• Charging: A bolt symbol next to the battery</li> </ul>	
9	<p>Chart tile</p> <p>The Chart tile offers a Numerical view and a Chart view.</p> 	<p>If Numerical view is visible:</p> <ul style="list-style-type: none"> <li>• Tap: Change to Chart view</li> </ul> <p>If Chart view is visible:</p> <ul style="list-style-type: none"> <li>• Tap: Change to Numerical view</li> <li>• Long press: Access slider to adjust the time axis</li> </ul>
10	<p>Statistic screen</p> <p>Shows the measurement statistics</p> 	<p>Swipe left / right to switch between Measurement and Statistic screen</p>

6


Note: Statistic Screen is only available on FlowMeter F2.

### 6.4.2 Configuration of the Measurement screen

The Measurement screen offers 1 chart tile and 3 value tiles which can be configured individually.



Long press the tile for 1 sec. to open the Edit parameter screen. Here you can select between all available values and units. Use the swipe up / down gesture to scroll through the lists of parameters.



All available values and units are listed in [→11 Measurement parameters and units](#)


### 6.4.3 Configuration of the Statistic screen

For each row, long press for 1 sec. to select any measurement value.

	Min	Max	Mean	
Flow	30.5	-50.0	130	55.0 l/min
Volume	400	13.3	550	234 ml
Pabs	12.5	-50.0	130	55.0 kPa
Pdiff	230	-50.0	130	55.0 mbar
Temper	21.3	-50.0	130	55.0 °C
Humidity	30.5	-50.0	130	55.0 %

6


Long press on a row for 1 sec. to open the Edit parameter screen. Use the swipe up / down gesture to scroll through the lists of parameters.



All available values and units are listed in [→11 Measurement parameters and units](#)

## 6.5 Menu


By tapping the **Menu** button, the sub menus are displayed:



No.		Description	Action
1	Zero Calibration	Opens the sub menu zero calibration ( <a href="#">→6.5.8 Zero offset calibration</a> )	Tap: The sub menu opens
2	Start Recording	<p>Opens the sub menu with the recording settings:</p> <ul style="list-style-type: none"> <li>• USB connection</li> <li>• Sampling interval</li> <li>• Recording duration</li> </ul> <p>Once the settings are made, the recording can be started by pressing the button "Start recording".</p> <p>The data is stored as comma separated value (CSV) file with ASCII encoding (Windows-1252 / ISO-8859-1 code page).</p> <p>Note: When a recording is active, this menu point changes its appearance as follows and opens the sub menu to stop the recording:</p> <div style="background-color: #cccccc; padding: 5px; border: 1px solid #ccc;"> <p><span style="color: red;">●</span> Recording... Tab to stop recoding</p> </div>	<p>If recording is inactive:</p> <ul style="list-style-type: none"> <li>• Tap: Sub menu of recording settings opens</li> </ul> <p>If recording is active:</p> <ul style="list-style-type: none"> <li>• Tap: Sub menu for stop recording opens</li> </ul>
3	Screenshot!	Exports a screenshot as Bitmap file (BMP) to a connected USB drive. Screenshot is only possible when a USB drive is connected.	Tap: Perform the Screenshot
4	Settings	<p>Opens the sub menu device settings:</p> <ul style="list-style-type: none"> <li>• Measurement</li> <li>• Trigger</li> <li>• Device</li> <li>• Bluetooth</li> <li>• About</li> </ul> <p>(<a href="#">→6.5.1 Settings</a>)</p>	Tap: The sub menu opens
5	Close Menu	Hides the overlaying Menu	Tap in the background: Hide the Menu

### 6.5.1 Settings

In the **Settings** menu, all device settings can be made:




6

No.		Description	Action
1	Back	Go back to the main Menu	
2	Measurement	Opens the sub menu for the measurement settings: <ul style="list-style-type: none"> <li>• Gas Standard</li> <li>• Gas Type</li> <li>• Time Axis</li> <li>• Filter Type</li> </ul> (→6.5.2 Measurement settings)	Tap: The sub menu opens
3	Trigger	Opens the sub menu for the trigger settings: <ul style="list-style-type: none"> <li>• Start/Stop Signal</li> <li>• Start/Stop Flow Threshold</li> <li>• Start/Stop Edge</li> </ul> (→6.5.4 Trigger settings)	Tap: The sub menu opens
4	Device	Opens the sub menu for the device settings: <ul style="list-style-type: none"> <li>• Screen Brightness</li> <li>• Screen Timeout</li> <li>• Language</li> <li>• Serial Interface</li> <li>• Factory Reset</li> </ul> (→6.5.5 Device settings)	Tap: The sub menu opens
5	Bluetooth	Opens the sub menu for the bluetooth configuration (→6.5.6 Bluetooth settings)	Tap: The sub menu opens
6	About	Opens the about screen with the following information: <ul style="list-style-type: none"> <li>• Serial Number</li> <li>• Software Version</li> <li>• Calibration Interval</li> <li>• Last Device Message</li> <li>• Software Update</li> </ul> (→6.5.7 Software update)	Tap: The about screen opens

### 6.5.2 Measurement settings

The **Measurement** menu contains all measurement related settings:



No.		Description	Action
1	Back	Go back to <b>Settings</b> menu	
2	Gas Standard	Opens the gas standards settings with the following available choices: <ul style="list-style-type: none"> <li>• STP</li> <li>• ATP</li> <li>• BTPS-A</li> <li>• CTP</li> </ul>	Tap: Selection of all gas standards opens
3	Gas Type	Opens the gas type setting with the following available choices: <ul style="list-style-type: none"> <li>• Air</li> <li>• Air/O<sub>2</sub>, O<sub>2</sub> concentration, can be selected</li> <li>• CO<sub>2</sub><sup>1</sup></li> <li>• N<sub>2</sub></li> </ul>	Tap: Selection of all gas types opens
4	Time axis	Opens a slider to adjust the time axis with the following available ranges: 2 sec., 5 sec., 10 sec., 15 sec., 30 sec., 60 sec., 120 sec.	Tap: Slider to adjust the time axis appears  If the slider is opened: <ul style="list-style-type: none"> <li>• Slide to the right: time range increases</li> <li>• Slide to the left: time range decreases</li> </ul>
5	Filter type	Opens the filter type setting with the following available choices: <ul style="list-style-type: none"> <li>• None</li> <li>• Low</li> <li>• Medium</li> <li>• High</li> </ul> (→6.5.3 Filter type)	Tap: Selection of all filter types opens



The selected gas type must match the measured gas. Improper setting of gas type or standard can lead to measurement errors of up to 60%.

<sup>1</sup> Only available for the model FlowMeter F2

### 6.5.3 Filter types

Numeric values for the graph are displayed when the screen is refreshed every 50 ms, but measurement takes place every 1 ms.

The sampling interval for measurement parameters is 1 ms. To reduce fluctuation of the reading and to make measurements easier to read, a filter can be applied..

The following four options are available:

- None (indication of the last reading without averaging)
- Low (average over 1 sec.)
- Medium (average over 2 sec.)
- High (average over 5 sec.)


Medium filter is the default setting.




This filtering of readings only impacts the values shown on the Flow-Meter display. During recording, only the raw and unfiltered readings are displayed.

### 6.5.4 Trigger settings

The **Trigger** menu contains all settings to configure the trigger (for more information about triggers see (→7 Measuring volume):




No.		Description	Action
1	Back	Go back to <b>Settings</b> menu	
2	Start Signal	Opens the start signal setting with the following available choices: <ul style="list-style-type: none"> <li>• Flow</li> <li>• Absolute Pressure (in Flow)</li> <li>• Differential Pressure</li> </ul>	Tap: Selection of start signal settings opens
3	Start Flow Threshold	Opens the start flow threshold configuration mask  The unit of the corresponding measurement value is used and can not be changed here. To change the unit, go back to the measurement screen and configure the measurement tile accordingly (→6.4.2 Configuration of the Measurement screen).	Tap: Mask to increase (+) and decrease (-) the start flow opens
4	Start Edge	Displays the current setting of the start edge	Tap on toggle button: toggle between “Rising Edge” and “Falling Edge”
5	Stop Signal	Opens the start signal setting with the following available choices: <ul style="list-style-type: none"> <li>• Flow</li> <li>• Absolute Pressure (in Flow)</li> <li>• Differential Pressure</li> </ul>	Tap: Selection of stop signal settings opens
6	Stop Flow Threshold	Opens the stop flow threshold configuration mask  The unit can not be changed here. For changing the unit go back to main menu and tap on the grey field (→6.4.2 Configuration of the Measurement screen).	Tap: Mask to increase (+) and decrease (-) the stop flow opens
7	Stop Edge	Displays the current setting of the stop edge	Tap on toggle button: toggle between “Rising Edge” and “Falling Edge”
8	Reset	Opens a window to perform the trigger factory reset (insert picture) 	Tap on Reset!: Device trigger configuration will be reset to factory defaults




### 6.5.5 Device settings

The **Device** menu contains all settings to configure the device:



6

No.		Description	Action
1	Back	Go back to <b>Settings</b> menu	
2	Screen Brightness	Opens a slider to adjust the brightness of the screen with the following available choices: 10 %, 20 %, 30 %, 40 %, 50 %, 60 %, 70 %, 80 %, 90 %, 100 %	Tap: Slider to adjust the screen brightness appears  If the slider is opened: <ul style="list-style-type: none"> <li>• Slide to the right: screen brightness increases</li> <li>• Slide to the left: screen brightness decreases</li> </ul>
3	Screen Timeout	Opens the screen timeout setting with the following choices: None, 1 minute, 5 minutes, 15 minutes, 30 minutes  When the time since the last user input is longer than the screen timeout setting, the screen brightness is reduced to a low level to increase the battery runtime. Any tap on the screen restores the regular brightness level.	Tap: Selection of all screen timeout settings opens
4	Language	Opens the language setting with the following available choices: English, German, Spanish, French	Tap: Selection of all languages opens
5	Serial interface (RS-232)	Opens the serial interface (RS-232) setting with the following choices: IMT protocol, IMT fast protocol, IMT express protocol, TSI 5300 protocol	Tap: Selection of the serial interfaces (RS-232) opens
6	Factory Default	Opens a window to perform the factory reset	Tap on <b>Reset!</b> : Device settings and screen configurations will be reset to factory defaults




### 6.5.6 Bluetooth settings


The FlowMeter App called IMT Analytics can be downloaded from the IMT Analytics website or from Google Play store or App store on your mobile device. Be sure, that Bluetooth is enabled.



In **Settings** menu on your FlowMeter select the sub menu **Bluetooth** to enable it. Pair your mobile device with your FlowMeter by entering the pairing code provided by the FlowMeter.




Choose the desired FlowMeter on the FlowMeter App to track the measurement values.



If connection is successful established the Bluetooth icon shown on your FlowMeter changed from Bluetooth enabled to Bluetooth connected:



In the FlowMeter App, a long press on a tile for 1 sec. opens the Edit parameter screen. Use the drop-down menu to select parameters and units.




### 6.5.7 Software update

The software update file can be downloaded from the IMT Analytics website. A USB flash drive with FAT32 formatting is required to install new software.


To update the FlowMeter software:

1. Copy the software update file to the root directory of a FAT32-formatted USB flash drive
2. Plug the USB flash drive into the USB-A host port
3. Switch on the FlowMeter
4. Tap the **Menu** button to open the main menu
5. Select **Settings**
6. Select **About**
7. Select **SW update**

The device verifies if a valid software update file is available on the connected USB flash drive. A progress bar is shown on the screen as follows:



If the verification was successful, the following message is shown on the screen:







If the verification was unsuccessful, an error message is shown on the screen. See details about possible error messages below.

8. Tap **Press to start the update!**


The software update starts. A progress bar is shown on the screen as follows:



Once the software update is loaded, the device reboots automatically. This is indicated on the screen as follows:



During the reboot, the screen turns white for a couple of seconds. After that, the following message is displayed on the screen:








Once the update is completed, the measurement screen is displayed

9. Disconnect the USB flash drive from the FlowMeter

The device is ready for use

If the verification of a valid software update file was unsuccessful, one of the following error messages is shown on the screen:

Error message displayed	Troubleshooting
<p>SW update Please insert USB device!</p> 	<ol style="list-style-type: none"> <li>1. Tap <b>Retry!</b> to try it again</li> <li>2. Ensure that a USB flash drive is connected to the USB-A host port (<a href="#">→5.4 Electrical interfaces</a>)</li> <li>3. Ensure the battery charge is &gt; 50%</li> <li>4. Format the USB flash drive using the file system <b>FAT32</b>, then re-download the software update file from the IMT Analytics website, save the software update file named Flow-MeterUpdate.bin in the root directory of the USB flash drive, disconnect the USB flash drive from your computer safely (Eject) and try to install the software again</li> <li>5. Repeat the process with a different USB flash drive</li> </ol>
<p>SW update Invalid file system on USB device!</p> 	<ol style="list-style-type: none"> <li>1. Tap <b>Retry!</b> to try it again</li> <li>2. Format the USB flash drive using the file system FAT32, then re-download the software update file from the IMT Analytics website, save the software update file named Flow-MeterUpdate.bin in the root directory of the USB flash drive, disconnect the USB flash drive from your computer safely (Eject) and try to install the software again</li> <li>3. Repeat the process with a different USB flash drive</li> </ol>
<p>SW update No update found on USB device!</p> 	<ol style="list-style-type: none"> <li>1. Tap <b>Retry!</b> to try it again</li> <li>2. Ensure that a valid software update file named FlowMeterUpdate.bin is in the root directory of the USB flash drive</li> </ol>
<p>SW update Update file is corrupted!</p> 	<ol style="list-style-type: none"> <li>1. Tap <b>Retry!</b> to try it again</li> <li>2. Ensure the USB flash drive is connected to the device while the software update verification is in progress</li> <li>3. Format the USB flash drive using the file system <b>FAT32</b>, then re-download the software update file from the IMT Analytics website, save the software update file named Flow-MeterUpdate.bin in the root directory of the USB flash drive, disconnect the USB flash drive from your computer safely (Eject) and try to install the software again</li> </ol>

Error message displayed	Troubleshooting
	<ol style="list-style-type: none"><li>1. Tap <b>Retry!</b> to try it again</li><li>2. Ensure the USB flash drive is connected to the device while the software update is in progress</li><li>3. Format the USB flash drive using the file system <b>FAT32</b>, then re-download the software update file from the IMT Analytics website, save the software update file named Flow-MeterUpdate.bin in the root directory of the USB flash drive, disconnect the USB flash drive from your computer safely (Eject) and try to install the software again</li></ol>



During installation make sure the battery charge is at least 50% or the device is plugged into the power supply with mains power connected.




While installation is running, the device must not be switched off!

### 6.5.8 Zero offset calibration

Follow the procedure described below to achieve the most accurate readings:


- Perform the zero offset calibration when the FlowMeter has warmed up. Warm up takes approximately 30 minutes.
- To perform zero offset calibration you must remove all connected tubes and protection caps from the device
- Tap on menu icon, select the sub menu **Zero Calibration**, and select: **Start Calibration**



6


#### Calibration is running

The blue bar increases from left to right in 5 sec.



#### Results


Calibration is successfully completed:



The device is now ready to use  
or

#### Calibration failed

The device is **not** ready to use.



If the calibration failed, ensure that all tubes and the protection caps are removed, and that there is no pressure or flow in the channel during zero offset calibration. In environments with high turbulence, close one of the flow ports using a protection cap.




Be aware that the correct calibration of a medical device is an important process. Incorrect calibration may lead to patient harm and / or increase the treatment duration. In case the measurement values seem implausible, check the measurement setup and measurement equipment. When in doubt about the performance of the FlowMeter, send it back to the manufacturer for calibration and do not change any settings on the medical device.



## 7 Measuring volume

### 7.1 General principle


The FlowMeter measures volumes by integrating the flow between the start and stop trigger events. Flow, Absolute Pressure (in Flow) and Differential Pressure can be used as trigger source.



### 7.2 Finding the correct trigger values

If you set a trigger for the first time, it is important to know the characteristic of the signal that will be used for the trigger (flow or pressure). It is therefore advisable to analyze that characteristic by using the chart view on the device or record the data to a USB flash drive first. By analyzing the graph afterwards, the appropriate trigger settings can be determined.

The following example shows a flow curve with a potential problem:



The first two circles indicate the triggers that should be used to measure the volume of this event. However, after the event a small false signal is visible that can be caused for example by switching of valves. In this case a pressure trigger should be used.

## 8 Servicing and care

### 8.1 Guidelines for servicing and care

Careful service in compliance with the instructions is essential for ensuring that FlowMeter operates safely and efficiently. Only components recommended by the manufacturer may be used.



It is absolutely essential to comply with the guidelines and service instructions issued by the various manufacturers.

### 8.2 Notes about changing parts



The service routines listed below may only be performed by persons who are familiar with FlowMeter. All further repair work may only be performed by authorised trained professionals. Please also observe the information issued by the various manufacturers.

### 8.3 Preventive cleaning and servicing routines


To ensure that your device operates with precision and reliability for as long as possible, it is essential to perform the following servicing routines regularly:

#### Every four weeks

For best accuracy, visually check the screens inside the flow ports. Inspect the screens for contamination by looking through each end of the flow channel while pointing the other end of the flow channel at a light source.

If the screen is contaminated, clean it according to the following procedure:

- Unscrew the flow port retention nut **1**
- Remove the flow port connector **2** and the screen **3**
- Wash the screen under running water and dry it
- Reinsert the screen and flow port into the device
- Install and tighten the flow port retention nut





If the screen can not be cleaned under running water, exchange it with a new one (→9.2 Parts).

Note: If the flow channel was observed to be clean in the past, contamination can easily be recognized.

**Every 12 months**

Factory calibration ensures reliable measurement. To calibrate the FlowMeter by the manufacturer, visit the website [www.easy-cal.com](http://www.easy-cal.com)



Only use the spare parts recommended by the manufacturer (→9 Accessories and spare parts).

**8.4 Cleaning**

For cleaning, only use agents which do not affect the functioning of the FlowMeter:

- Flow channel: dry, compressed and clean air
- Enclosure: soft, lint-free cloth, dry or moist (soap, mild cleaning agent or iso-propyl alcohol)

For cleaning, do not use any solvent or chlorine based or abrasive detergents. Solvents can permanently offset the humidity sensor.

The device is not intended to be disinfected or reprocessed and shall therefore be kept clean.

Using incorrect cloths can cause scratches on the device.

**8.5 Trouble shooting**

Fault	Possible cause	Trouble shooting
FlowMeter not turning on	No power	Connect power supply, check voltage, check all connections, use a compliant power supply as described in the chapter specification
Error message: Sensor Failure!	Sensor defective or extreme environment condition Error code: xxxx	Follow allowed operation conditions and restart the device. If the error persists, contact the distributor or manufacturer for repair of the device
Error message: Battery not available!	Battery disconnected or defective	Restart the device and charge the battery for at least 2 hours. If the error persists, contact the distributor or manufacturer for repair of the device
Error message: Measurement not available	Measurement subsystem not responding	Restart the device. If the error persists, contact the distributor or manufacturer for repair of the device
Error message: Calibration overdue!	12 month since last calibration	Contact IMT Analytics for calibration: <a href="http://www.easy-cal.com">www.easy-cal.com</a>

Fault	Possible cause	Trouble shooting
Error message: Calibration is due	11 month since last calibration	Contact IMT Analytics soon for calibration: <a href="http://www.easy-cal.com">www.easy-cal.com</a>
Error message: Device over-heated	Channel temperature > 65 °C or Battery temperature > 60 °C or any other temperature > 75 °C	Move device to cooler environment to avoid damage.  Any further increase in temperature will damage the device!
Error message: Channel pressure too high	Abs. Channel pressure > 2.05 bar	Reduce channel pressure below 2 bar absolute pressure to avoid damage.  Any further increase in pressure will damage the device!
Error message: Ambient pressure too high	Abs. ambient pressure > 1250 mbar	Reduce ambient pressure below 1200 mbar absolute pressure to avoid damage.  Any further increase in pressure will damage the device!
Error message: Ambient pressure too low	Abs. ambient pressure < 490 mbar	Increase the ambient pressure above 540 mbar absolute pressure to avoid damage.  Any further decrease in pressure will damage the device!
Error message: External differential pressure too high	Ext DP sensor < -250 mbar or Ext DP sensor > 250 mbar	Reduce differential pressure to ±250 mbar to avoid damage.  Any further increase in pressure will damage the device!
The display values of the FlowMeter do not match the device under test (DUT)	<ul style="list-style-type: none"> <li>Flow channel dirty</li> <li>Calibration due</li> <li>Error in setup measurement</li> </ul>	<ul style="list-style-type: none"> <li>Visually inspect the flow channel</li> <li>Contact IMT Analytics soon for calibration: <a href="http://www.easy-cal.com">www.easy-cal.com</a></li> <li>Check measurement setup</li> </ul>
Touch screen not responsive	Software locked up	Hard reset the device. To do so: <ol style="list-style-type: none"> <li>1. Disconnect any cables including the power supply</li> <li>2. Press and hold the power button for at least 10 sec.</li> <li>3. Device will shut down and restart automatically</li> </ol>

## 8.6 Contact

If you have any questions or problems please contact one of the departments listed below.

### 8.6.1 Manufacturer's name and address

IMT Analytics AG  
Gewerbstrasse 8  
CH-9470 Buchs  
Switzerland

Tel: +41 (0)81 750 67 10

Email: [sales@imtanalytics.com](mailto:sales@imtanalytics.com)

### 8.6.2 Technical support

Tel: +41 (0)81 750 67 10

Email: [techsupport@imtanalytics.com](mailto:techsupport@imtanalytics.com)

## 9 Accessories and spare parts



### 9.1 Ordering address

IMT Analytics AG  
 Gewerbestrasse 8  
 CH-9470 Buchs  
 Switzerland

Tel: +41 (0)81 750 67 10  
 Email: [sales@imtanalytics.com](mailto:sales@imtanalytics.com)

### 9.2 Parts

Picture	Name	Article Number
	Power Supply including mains adapters for different countries	700.421.000
	Protection Filter RT019	302.531.000
	Flow Port Screen	700.412.000
	Adapter set	300.548.000
	RS-232 Interface Cable	302.075.000

Picture	Name	Article Number
 <p>A grey carrying bag with a handle, featuring the text 'IMT.Analytics' and 'FlowMeter' on its front panel.</p>	<p>Carrying Bag FlowMeter</p>	<p>700.422.000</p>
 <p>A blue component with a white cap and a metal ring, used for protecting the flow channel. The blue part has markings '0' and '1'.</p>	<p>Flow Channel Protection Cap</p>	<p>302.780.000</p>

## 10 Disposal

Disposal of the device is the operator's responsibility. Keep in mind the FlowMeter contains a Li-Ion battery and must be disposed of accordingly.

The device can:

- be delivered, carriage free and duty paid to the manufacturer for disposal.
- be handed over to a licensed private or public collection company.
- be professionally broken down into its constituent parts by the operator and be recycled or disposed of in accordance with regulations.

In the case of self-disposal the disposal regulations are country-specific and are contained in relevant laws and ordinances. These codes of conduct must be obtained from the authorities responsible.

In this context, waste material must be recycled or destroyed

- without endangering human health
- without using processes or methods that harm the environment, especially water, air, soil, animals and plants
- without causing noise or odour nuisances
- without having a detrimental effect on the environment or landscape.



## 11 Measurement parameters and units

### 11.1 Measurement parameters and units

#### 11.1.1 Pressure readings

Measurement parameter	Designation	Units of measurement
Absolute pressure in flow channel	P <sub>Abs</sub>	mbar, inH <sub>2</sub> O, cmH <sub>2</sub> O, PSI, mmHg, hPa, kPa, Pa
Differential pressure	P <sub>Diff</sub>	

#### 11.1.2 Flow and Volume readings

Measurement parameter	Designation	Units of measurement
Flow	Flow	L/min, ft <sup>3</sup> /min
Volume	Volume	L, ml, ft <sup>3</sup>

#### 11.1.3 Temperature and Humidity readings

Measurement parameter	Designation	Units of measurement
Temperature	Temperature	°C, K, °F
Humidity	Humidity	% RH

## 12 Wireless communication

FlowMeter contains a generation 5.1 bluetooth module for wireless communication. The module complies with all requirements from EN 300 328 V2.2.2.

FCC ID: QQQGM210P  
 IC ID: 5123A-GM210P

The module is restricted to the following channels and output power.

Band	Lowest frequency	Highest frequency	Maximum power (EIRP) in dBm
2.4 GHz	2400 MHz	2483.5 MHz	10

FCC RF Radiation Exposure Statement:

This the equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment.

End users must follow specific operating instructions for satisfying RF exposure compliance:

- The transmitter meets the Mobile requirements at a distance of 20 cm and above from the human body, in accordance with the limit(s) exposed in the RF Exposure Analysis
- The transmitter also meets the Portable requirements at distances equal or above 5.3 mm\*

Part 15 information to the user:

The equipment has been tested and found to comply with the limits of Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and the receiver
- Connect the equipment into a outlet on a circuit different from that to which the receiver is connected
- Consult dealer or an experienced TV technician for help

NOTE: "Harmful interference" is defined in 47 CFR §2.122 by the FCC as follows: Interference which endangers the functioning of a radionavigation service or of other safety services or seriously degrades, obstructs, or repeatedly interrupts a radio communication service operating in accordance with [the ITU] Radio Regulations.

\* The mechanical design of the FlowMeter ensures the distance to the transmitter is always equal or above 5.3 mm

## 13 Appendix

### 13.1 Abbreviations and glossary

<b>A</b>	
A	Ampere
<b>B</b>	
bar	1 bar = 14.50 psi
<b>C</b>	
°C	Degrees Celsius Conversion from Celsius (C) to Fahrenheit (F): $F = 9 \times C / 5 + 32$
cm	Centimeters
cmH <sub>2</sub> O	Centimeters of water column
<b>D</b>	
dBm	Decibel-milliwatts
<b>F</b>	
°F	Degrees Fahrenheit Conversion from Fahrenheit (F) to Celsius (C): $C = (F - 32) \times 5 / 9$
ft <sup>3</sup>	Cubic foot
ft <sup>3</sup> /min	Cubic foot per minute
<b>G</b>	
GHz	Gigahertz
GND	Ground
<b>H</b>	
Hz	Hertz (1 Hz = 1 / sec.)
h	Hour
hPa	Hectopascal
<b>I</b>	
IP	Protection class according to standard
inH <sub>2</sub> O	Inches of water column
<b>K</b>	
K	Kelvin
kPa	Kilopascal
<b>L</b>	
L	Litre
L/min	Litres per minutes
<b>M</b>	
MB	Megabyte
mbar	Millibar
MHz	Megahertz
mL	Millilitre
mm	Millimetre
mmHg	Millimeters of Mercury
ms	Millisecond

**P**

Pa	Pascal
psi	Pounds per square inch (1 bar = 14.50 psi)

**R**

RH	Relative Humidity
RS-232	Serial interface
RJ-10 FCC	Connector for serial interface (telephone connector according to FCC registration, U.S. Federal Communications Commission; RJ = 'Registered Jack')

**S**

sec.	seconds
------	---------

**V**

$\dot{V}$	Flow
V	Volt
VAC	Volts Alternating Current
VDC	Volts Direct Current

## 13.2 Index

**A**

Abbreviations and glossary 52  
Accessories 46  
Acoustic feedback 22  
Appendix 52  
Approvals 13

**B**

Battery operation 13  
Bluetooth settings 34

**C**

Calibration by user 11  
Changing parts 42  
Cleaning 43  
Communication interfaces 10  
Contact 43, 45

**D**

Density 10  
Device labels 14, 15  
Differential pressure 19  
Directives 13  
Disposal 48, 49, 50  
Dynamic viscosity 10

**E**

Electrical interfaces 20, 37

**F**

Feature Set 15  
Filter types 31  
Finding the correct trigger values 41  
Flow and Volume readings 49  
Flow channel 18  
Flow readings 49

**G**

Gas standards 12  
General principle 41  
Gestures 22  
Glossary 52  
Guarantee 7

**I**

Intended use 6  
Introduction 5

**M**

Measurement 24  
Measurement parameters 9  
Measurement parameters and units 49  
Measurement screen 24  
Measurement settings 30  
Measuring volume 41  
Mechanical connectors 17  
Menu 28

**O**

Operating data 11  
Operation 22

**P**

Personnel 7  
Physical data 10  
Power supply 12, 17  
Pressure readings 49  
Preventive cleaning 42  
Protection Filter RT019 17

**R**

Responsibility 7  
Responsibility and guarantee 7  
RS-232 20

**S**

Safety instructions 6, 7  
Screens 24  
Service life 8  
Servicing 42  
Servicing and care 42  
Servicing routines 42  
Settings 28, 29  
Software updates 34, 35  
Spare parts 46  
Special functions 10  
Specifications 9  
Start-up 16  
Supply voltage 17  
Switching the device on and off 22  
Symbols 14

**T**

Trigger settings 29, 32  
Trouble shooting 43

**U**

USB 21  
USB to Computer 20  
User control 22, 24  
User control gestures 22, 24

**V**

Version information 5

**W**

Wireless communication 50

**Z**

Zero Calibration 28, 39

# IMT.Analytics

IMT Analytics AG . Gewerbestrasse 8 . 9470 Buchs . Switzerland  
T +41 81 750 67 10 . [www.imtanalytics.com](http://www.imtanalytics.com)