

SmartLung 2000

Test Lungs adjustable

testlung

the unique simulation

SmartLung 2000 2L – the cleverest and most reliable way to check function and precision of ventilators and anaesthesia machines

Offers settings according to standards IEC 60601-2-12, IEC 60601-2-13, EN 794-3

- Simulates adult lungs
- The only test lung on the market which covers all the necessary checkpoints
 - Compliance, Resistance and leakage widely adjustable
 - Compliance: 200, 500 mL/kPa
 - Compliance: 25, 60, 75 mL/mbar
 - Resistance 5, 20, 50, 200 mbar/L/s
- Extremely handy and user-friendly
- Top performance and quality at a low price
- Incl. Carrying Bag



Features SmartLung 2000

Easy to Handle

Practical transport bag.

Adjustable Leakage

Leakage can be adjusted by turning the side screw. No adapter required.

Adjustable Resistance

Different airway resistances can be simulated simply by turning the connector.

Adjustable Compliance

The slider enables different degrees of lung compliance to be simulated without an adapter.



IMT.Analytics

IMT Analytics AG . Gewerbestrasse 8 . 9470 Buchs . Switzerland
T +41 81 750 67 10 . www.imtanalytics.com